





---

ACBS ☒ GMS ☒

---

**Description****Food for special medical purposes.**

MCTprocal is a neutral tasting powder high in medium chain triglycerides (MCT).

**Indications**

MCTprocal is for the dietary management of disorders of long-chain fatty acid oxidation, fat malabsorption and other disorders requiring a high MCT, low long-chain triglyceride (LCT) diet from 3 years onwards.

**Dosage and Administration**

To be determined by the clinician or dietitian and is dependent on the age, body weight and medical condition of the patient.

**MCTprocal (16g) = 10g MCT, 112kcal and 2g protein.**

**Preparation Guidelines**

Add MCTprocal to the desired food or drink, stir until the powder is dissolved.

MCTprocal can also be used in cooking and baking.

Not suitable for addition to fruit juice.

Serving suggestions for MCTprocal are available on request.

**IMPORTANT NOTICE**

Use under medical supervision.

Not suitable for use as a sole source of nutrition.

Suitable from 3 years of age and onwards.

For enteral use only.

## Nutritional Information

		per 100g	per 16g sachet			per 100g	per 16g sachet
Energy	kJ	2907	465	<b>Minerals</b>			
	kcal	703	112				
Fat	g	63.5	10	Sodium	mg	220	35
	of which saturates	g	61.2		mmol	9.5	1.5
MCT	g	60.7	10	Potassium	mg	610	98
					mmol	15.3	2
Carbohydrate	g	20.6	3.3	Calcium	mg	360	58
	of which sugars	g	3.1		mmol	9	1.4
Protein	g	12.2	2.0	Phosphorus	mg	500	80
Salt	g	0.55	0.09		mmol	16	2.6
				<b>Others</b>			
				LCT	g	0.90	0.1

### Osmolality:

1 sachet + 50ml water = 198mOsm/kg

### Ingredients:

Medium Chain Triglyceride (MCT) oil, Dried Glucose Syrup, Sodium Caseinate (**Milk**), Stabilizer (E340), Anti-caking Agent (E341)

### Allergen Declaration

Contains **Milk** (milk protein).

### Storage

**Unopened:** Store in a cool, dry place.

**Opened:** Sachets are designed for single use and should be used once opened. Any unused powder in an open sachet should be kept in an airtight container or sealed with a clipper and used within 24 hours.

### Pack Size/Weight

30 x 16g e sachets = 480g

### Sample pack

12 x 16g e sachets = 192g



Enhancing Lives Together